

<u>TEST NAME</u>	<u>TEST COD E</u>	<u>METHOD</u>	<u>SPECIME N</u>	<u>QTY</u>	<u>SAMPLE COLLECT ION</u>	<u>REPORT</u>
ALANINE AMONOTRANSFERASE (SGPT)	S237	spectrophotometry	serum	2 ml	Red top	2 hours
ALBUMIN	S107	spectrophotometry	serum	2 ml	Red top	2 hours
ALBUMIN+GLOBULIN+ A/G RATIO	S315	spectrophotometry	serum	2 ml	Red top	2 hours
ALKALINE PHOSPHATASE	S108	spectrophotometry	serum	2 ml	Red top	2 hours
AMMONIA	S110	spectrophotometry	WHOLE BLOOD EDTA	2 ml	PURPLE top	2 hours
AMYLASE,24HRS URINE	S316	spectrophotometry	24 hours urine	-	Urine container	2 hours
AMYLASE serum	S111	spectrophotometry	serum	2 ml	Red top	2 hours
ANION GAP	S310	Spectrophotometry, ISE	serum	2 ml	Red top	2 hours
ASPARATE AMINOTRANSFERASE (SGOT)	S236	spectrophotometry	serum	2 ml	Red top	2 hours
BILIRUBIN (T,D,I)	S400	spectrophotometry	serum	2 ml	Red top	2 hours
BILIRUBIN, DIRECT	S124	spectrophotometry	serum	2 ml	Red top	2 hours
BILIRUBIN, TOTAL	S123	spectrophotometry	serum	2 ml	Red top	2 hours
BLOOD UREA NITROGEN	S136	spectrophotometry	serum	2 ml	Red top	2 hours
CALCIUM,24HRS URINE	S273	spectrophotometry	24 hours urine	2 ml	Urine container	2 hours
CALCIUM total SERUM	S139	spectrophotometry	serum	2 ml	Red top	2 hours
CALCIUM/CREATININE RATIO, URINE	S274	spectrophotometry	urine	2 ml	Urine container	2 hours
CHLORIDE SERUM	S145	spectrophotometry	serum	2 ml	Red top	2hours
CHLORIDE urine	S275	spectrophotometry	urine	2 ml	Urine container	2hours
CHOLESTEROL HDL	S187	spectrophotometry	serum	2 ml	Red top	2 hours
CHOLESTEROL TOTAL	S146	spectrophotometry	serum	2 ml	Red top	2 hours
CREATININE KINASE (CPK)	S151	spectrophotometry	serum	2 ml	Red top	2 hours

CK (MB)	S152	spectrophotometry	SERUM	2ml	Red Top	2 hours
ABG	S102	ISE	WB HEPARINI SED	1ML	Heparinis ed syringe	1hr
ADA ( Fluid / serum)	S105	spectrophotometry	Fluid	2ml	Red top	3hrs
Ammonia blood	S110	spectrophotometry	Whole EDTA blood	2ml	Purple top	2hrs
Blood Ketones	S132	Amperometric blood chemistry	Whole blood	1ml	-	1hr
Urea, Blood	S136	spectrophotometry	SERUM	2ml	Red Top	2 hours
CALCIUM SERUM TOTAL	S139	spectrophotometry	SERUM	2ml	Red Top	2 hours
Calcium ionized	S196	POTENTIOMETRY	Herapinis ed WB	2ml	Heparinis ed SYRINGE	1hour
C Reactive protein (quantitative)	S154	immunoturbidimetr y	SERUM	2ml	Red Top	2 hours
CREATININE CLEARANCE,	S150	spectrophotometry	SERUM	2 ml	Red top & Urine container	2 hours
CREATININE,24HRS URINE	S276	spectrophotometry	24 hours urine	24 hour urine	Urine Container	2 hours
CREATININE SERUM	S153	spectrophotometry	serum	2 ml	Red top	2 hours
DIRECT LDL CHOLESTEROL	S200	spectrophotometry	serum	2 ml	Red top	2 hours
ELECTROLYTES (NA/K)	S241	spectrophotometry	serum	2 ml	Red top	2 hours
Electrolytes total	S262	spectrophotometry	serum	2 ml	Red top	2 hours
GAMMA GLUTAMYL TRANSFERANCE	S179	spectrophotometry	serum	2 ml	Red top	2 hours
GLUCOSE FASTING	S133	spectrophotometry	serum	2 ml	Gray top	2 hours
GLUCOSE POST- PRANDIAL	S134	spectrophotometry	serum	2 ml	Gray top	2 hours
GLUCOSE RANDOM	S135	spectrophotometry	serum	2 ml	Gray top	2 hours
GLUCOSE TOLERANCE TEST-4 SAMPLE	S305	spectrophotometry	serum	2 ml	Gray top & urine container	2 hours

GLUCOSE TOLERANCE TEST (GTT-2)	S304	spectrophotometry	serum	2 ml	Gray top & urine container	2 hours
GLUCOSE, URINE	S292	Strip method/ manual	Random Urine	2 ml	Urine container	
GLYCOSYLATED HEMOGLOBIN (HbA1C)	S181	spectrophotometry	EDTA whole blood	2 ml	Purple Top	2 hours
IRON, S	SR3532	spectrophotometry	serum	2 ml	Red top	2 hours
LACTATE DEHYDROGENASE	S199	spectrophotometry	Serum/fl uid	2 ml	Red top, container	2 hours
MAGNESIUM SERUM	S205	spectrophotometry	serum	2 ml	Red top	2 hours
MICROALBUMINURIA, URINE		spectrophotometry	urine	2 ml	Sterile container	2 hours
GGT	S179	spectrophotometry	urine	2 ml	Sterile container	2 hours
OSMOLALITY, SERUM	S311	spectrophotometry	serum	2 ml	Red top	2 hours
LIPASE , SERUM	S204	spectrophotometry	urine	2 ml	Sterile container	2 hours
SODIUM/ POTASSIUM , URINE	S285	ISE	URINE	2ML	URINE CONTAINER	2HRS
PHOSPHORUS, 24HRS URINE	S289	spectrophotometry	urine	2 ml	Urine container	2 hours
PHOSPHORUS, S	S217	spectrophotometry	serum	2 ml	Red top	2 hours
POTASSIUM, S	S219	ISE	serum	2 ml	Red top	2 hours
POTASSIUM URINE	S283	ISE	serum	2 ml	Red top	2 hours
PROTEIN/CREATININE RATIO	S291	spectrophotometry	serum	2 ml	Red top	2 hours
PROTEIN, 24HRS URINE	S101	spectrophotometry	urine	2 ml	Sterile container	2 hours
SAAG	S235	spectrophotometry	Serum & Astic Fluid	2 ml	Red Top + Sterile container	2 hours
SODIUM, S	S239	ISE	Serum	2 ml	Red top	2 hours
SODIUM, URINE	S240	I SE	urine	2 ml	Urine Container	
SODIUM & POTASSIUM, serum	S241	I SE	Serum	2 ml	Red top	2 hours
TOTAL PROTEIN, FLUID	S222	spectrophotometry	Fluid	2 ml	Plain vial	2 hours
TOTAL PROTEIN, serum	S222	spectrophotometry	serum	2 ml	Red top	2 hours



## 24 HRS WORKING

663 GTB NAGAR, NEAR GURU RAVIDASS CHOWK  
 JALANDHAR (Phone Lab.0181 2463300,Recp.01814679999)

### DIRECTORY OF TESTS RATE w.e.f.- 01/04/2014

S.NO.	TEST NAME	Minimum sample requirement	Vial /container	Additive	METHOD	TAT	STABILITY/storage
1	Absolute Eosinophil count (AEC)	2.0 ml whole blood	Purple top tube	EDTA	Cell counter & microscopy	2 hours	Room temp. 6 hours and 48 hours at 4-8 degree Celcius
2	ABG	0.2ml arterial blood	Heparinised syringe	Heparin (1000IU/ml)	ISE	30 min	Immediately.
3	Adenosine Deaminase (ADA)	1.0 ML serum plasma, body fluid and CSF	Red top tube	Nil	Spectrophotometry	4 hours	Room temp. 2 hours and 07 days at 4-8 degree Celcius
4	AFB STAIN	Sputum, urine,Pus, body fluid FNAC &Histo Slides	Sterilize container	NA	Microscopy	01 day	Room temp. 2 hours , 48 hours at 4-8 degree Celcius
5	AFB CULTURE	Sputum, tissue in normal saline/ whole blood heparinised	Sterilize container	NA	LJ	01 Month	Room temp. 2 hours , 48 hours at 4-8 degree Celcius
6	Albumin, Serum	0.5 ml serum	Red top tube	Nil	Spectrophotometry	01 hour	Room temp. 2 hours , 07 days at 4-8 degree Celcius
7	Alkaline phosphatase, Serum	0.5 ml serum	Red top	Nil	Spectrophotometry	01 hour	Room temp. 72 hours and 07 days at 4-8 degree Celcius
8	Ammonia	2ml	Purple top tube	EDTA	Spectrophotometry	2 hours	immediately
9	Amylase	0.5 ml serum	Red top	Nil	Spectrophotometry	01 hour	Room temp. 6 hours and 07 days at 4-8 degree Celcius
10	Anti-StreptolysinASO	0.5 ml serum	Red top	Nil	Slide test	01 hour	Room temp. 4-6 hours and 02 days at 4-8 degree Celcius

TOTAL PROTEIN,ALBUMIN,GLOBALULIN		spectrophotometry	serum	2 ml	Red top	2 hours
TRIGLYCERIDES	S264	spectrophotometry	serum	2 ml	Red top	2 hours
UREA NITROGEN,24HRS URINE	S318	spectrophotometry	24 HOURS URINE		Urine container	2 hours
URIC ACID.24HRS URINE	S319	spectrophotometry	24 HOURS URINE		CONTAINER	2 hours
URIC ACID	S271	spectrophotometry	serum	2 ml	Red top	2 hours
antiTTG		ELISA	serum	2 ml	Red top	2 DAYS
LACTATE	S198	Potentiometry	Herapinised /EDTA WB	2ml	Green/ purple top	1hour
LDL ( DIRECT )	S200	spectrophotometry	serum	2 ml	Red top	2 hours
LIPASE,serum	S204	spectrophotometry	serum	2 ml	Red top	2 hours
TROPONIN T	S265	RAPID CARD TEST	WB/serum/plasma	2ml	Red /PURPLE top	1hrs
TROPONIN I	S264	RAPID CARD TEST	WB/serum/plasma	2ml	Red/ purple top	1hrs

## CLINICAL PATHOLOGY

<u>TEST NAME</u>	<u>TEST CODE</u>	<u>METHOD</u>	<u>SPECIMEN</u>	<u>MIN QNT</u>	<u>CONTAINER</u>	<u>TAT</u>
STOOL FOR ISOSPORA	S248	microscopy	stool	2g	STERILE CONTAINER	2hrs
STOOL FOR FAT GLOBULES	S246	microscopy	stool	2g	STERILE CONTAINER	2hrs
STOOL HANGING DROP	S247	microscopy	stool	2g	STERILE CONTAINER	2hrs
STOOL FOR REDUCING SUBSTANCE	S250	manual	stool	2g	STERILE CONTAINER	2hrs
Montoux test	S208	manual	-	-	-	48 hrs
SEMEN ANALYSIS	S233	PHYSICAL/ MICROSCOPY	SEMEN	1 ML	STERILE CONTAINER	01 HOURS
STOOL FOR OCCULT BLOOD	S212	Rapid card	STOOL	05 GM	STERILE CONTAINER	02 HOURS
STOOL ROUTINE examination	S244	MICROSCOPY	STOOL	05 GM	STERILE CONTAINER	02 HOURS
STOOL FOR CRYPTOSPORidium	S245	MICROSCOPY	STOOL	05 GM	STERILE CONTAINER	02 HOURS
URINALYSIS ROUTINE	S231	DIPSTICK/ MICROSCOPY	URINE	15 ,ML	STERILE CONTAINER	01 HOURS
URINE FOR BILE PIGMENTS	S268	MANUAL/ Distich	URINE	15 ,ML	STERILE CONTAINER	01 HOURS
URINE FOR BILE SALT	S267	MANUAL	URINE	15 ,ML	STERILE CONTAINER	01 HOURS
URINE FOR REDUCING SUBSTANCES	S284	MANUAL	URINE	15 ,ML	STERILE CONTAINER	01 HOURS
URINE FOR UROBILONOGEN	S287	MANUAL/ Distick	URINE	15 ,ML	STERILE CONTAINER	01 HOURS
URINE FOR FUNGAL HYPHAE	S279	microscopy	urine	5ml	STERILE CONTAINER	1hr
URINE FOR FAT GLOBULES	S278	microscopy	urine	5ml	STERILE CONTAINER	1hr
URINE FOR Hb	S280	MANUAL/ Distick	URINE	15 ,ML	STERILE CONTAINER	01 HOURS
URINE FOR NITRITE	S281	MANUAL/ Distick	URINE	15 ,ML	STERILE CONTAINER	01 HOURS
BJ PROTIEN	S117	rapid	urine	5ml	STERILE CONTAINER	1hr
URINE PORPHOBILINOGEN	S308	MANUAL	urine	5ml	STERILE CONTAINER	2hr

URINE KETONES	S288	MANUAL	urine	5ml	STERILE CONTAINER	2hr
URINE SPECIFIC GRAVITY	S286	MANUAL	urine	5ml	STERILE CONTAINER	2hr
URINE PH	S282	MANUAL	urine	5ml	STERILE CONTAINER	2hr
Stool for ph	S249	MANUAL	stool	5g	STERILE CONTAINER	1hr
Routine semen examination	S233	MANUAL	semen	5ml	STERILE CONTAINER	2hr
TESA	S256	MICROSCOPY	TESTICULAR TISSUE	5ml	STERILE CONTAINER	2hr

## HEMATOLOGY

<u>TEST NAME</u>	<u>TEST CODE</u>	<u>METHOD</u>	<u>SPECIMEN</u>	<u>QTY</u>	<u>SAMPLE COLLECTION</u>	<u>REPORT</u>
LE CELL	S201					
ABO GROUP & RH TYPE	S131	SLIDE AGGLUTINATION	EDTA BLOOD	2 ML	EDTA	01 HOUR
ABS EOSINOPHIL COUNT	S103	MICROSCOPY	EDTA BLOOD	2 ML	EDTA	01 HOUR
ABS NEUTROPHIL COUNT	S104	MICROSCOPY	EDTA BLOOD	2 ML	EDTA	01 HOUR
APTT	S113	Coagulometer	CITRATE BLOOD	2 ML	CITRATE TUBE	01 HOUR
BT & CT	S138	MANUAL	OTHERS	OTHERS	OTHERS	30 MINUTES
BONE MARROW EXAMINATION	S137	MICROSCOPY	BONE MARROW SMEAR/EDTA WHOLE BLOOD	SMEARS/3 ML in heparin	SLIDES/EDTA	NEXT DAY
COMPLETE BLOOD COUNT (CBC)	S141	AUTO ANALYSER	EDTA BLOOD	2 ML	EDTA	01 HOUR
ESR	S164	WESTERGEREN	CITRATE WHOLE BLOOD	2 ML	BLACK TOP	02 HOUR
HAEMOGLOBIN (HB)	S188	Autoanalyser/spectrophotometer	EDTA BLOOD	2 ML	EDTA	01 HOUR
HB,TLC & DLC		AUTO ANALYSER	EDTA BLOOD	2 ML	EDTA	01 HOUR
HB,TLC,DLC &PLT	S185	AUTO ANALYSER	EDTA BLOOD	2 ML	EDTA	01 HOUR

HEMOGLOBIN & HEMATOCRIT		AUTO ANALYSER	EDTA BLOOD	2 ML	EDTA	01 HOUR
OSMOTIC FRAGILITY	S213	MANUAL	HEPARIN WHOLE BLOOD	04 ML	GREEN TOP	24 hrs
SICKLING PHEN.	S238	Microscopy	EDTA BLOOD	2 ML	EDTA	2 hrs
G6PD	S178	manual	EDTA BLOOD	2 ML	EDTA	2 hrs
COOMBS TEST DIRECT	S148	manual	Edta wb	2ml	Purple top	3hrs



COOMBS TEST INDIRECT	S147	manual	serum	2ml	Red top	3hrs
MALARIAL Parasite smear	S207	MICROSC OPY	EDTA BLOOD	2 ML	EDTA	02 HOUR
PBF	S216	MICROSC OPY	EDTA BLOOD	2 ML	EDTA	02 HOUR
PLATELET COUNT	S218	MICROSC OPY+ autoanalys er	EDTA BLOOD	2 ML	EDTA	01 HOUR
PROTHROMBIN TIME ( PTI)	S224	COAGULO METER	Citrate Blood	2 ml	Blue Vail	2 Hrs
RETICULOCYTE COUNT	S228	MICROSC OPY	EDTA BLOOD	2 ML	EDTA	02 HOUR
SPUTUM /URINE EOSINOPHILS	S309	MICROSC OPY	SPUTUM/U RINE	2 ML	STERILE CONTAINER	012HOUR
WHITE CELL COUNT (TLC)	S261	autoanalys er	EDTA BLOOD	2 ML	EDTA	01 HOUR
Packed cell volume	S214	autoanalys er	EDTA BLOOD	2 ML	EDTA	01 HOUR
Mucin clot test	S209	manual	Synovial fluid	5ml	Sterile container	3hrs

## CYTOLOGY AND HISTOPATHOLOGY

<u>TEST NAME</u>	<u>TEST CODE</u>	<u>METHOD</u>	<u>SPECIMEN</u>	<u>QTY</u>	<u>SAMPLE COLLECTION</u>	<u>REPORT</u>
FLUID for malignant cells	S170	Microscopy	Any Body fluids	5 ml	Sterile container	Same day
FLUID for CYTOLOGY	S168	Microscopy	ANY Body fluids	5 ml	Sterile container	Same day
BAL FLUID FOR CYTOLOGY	S119	Microscopy	BAL fluid	5 ml	Sterile container	Same day
FNAC	S171	MANUAL	SMEARS/S LIDES	PREPARED SMEARS	SMEARS	SAME DAY
LARGE SPECIMEN	S190	MANUAL	TISSUE	WHOLE TISSUE	TISSUE IN 10 % Formalin with complete form	7 days
Small specimen	S191	MANUAL	TISSUE	WHOLE TISSUE	TISSUE IN 10 % Formalin with complete form	7 days
PAP smear	S215	manual	Smears/slides	Prepared slides	slides	Same day
Histo/ Cyto For review	S229/S 230	manual	Blocks & slides	Prepared slides & blocks	Slides/ Blocks	2 days

## HORMONES/VITAMINS/TUMOR MARKERS/ANEMIA MARKERS

<u>TEST NAME</u>	<u>TEST CODE</u>	<u>METHOD</u>	<u>SPECIMEN</u>	<u>QTY</u>	<u>SAMPLE COLLECTION</u>	<u>REPORT</u>
AMH	S109	ECLIA	SERUM	2ml	Red top	2 hours
AFP	S320	ECLIA	Serum	2ml	Red top	2 hours
ANTI TPO	S321	ECLIA	Serum	2ml	Red top	2 hours
BETA-HCG	S121	ECLIA	SERUM	2ml	Red top	2 hours
DHEAS	S160	ECLIA	Serum	2ml	Red top	2 Hours
Growth Hormone	S183	ECLIA	Serum	2ml	Red top	2 Hours
IgE	S323	ECLIA	SERUM	2ml	Red top	2 hours
ESTRADIOL	S165	ECLIA	Serum	2 ml	Red top	2 hours
FOLLICLE STIMULATING HORMONE (FSH)	S174	ECLIA	Serum	2 ml	Red top	2 hours
FREE THYROXINE (FT4)	S173	ECLIA	Serum	2 ml	Red top	2 hours
FREE TRIIODOTHYRONINE (FT3)	S172	ECLIA	Serum	2ml	Red top	2 hours
CA125	S308	ECLIA	Serum	2ml	Red top	2 hours
Insulin ,serum	S232	ECLIA	Serum	2ml	Red top	2 Hours
LUTEINIZING HORMONE (LH)	S202	ECLIA	Serum	2 ml	Red top	2 hours
PROLACTIN	S221	ECLIA	Serum	2 ml	Red top	2 hours
PROSTATE SPECIFIC ANTIGEN	S223	ECLIA	Serum	2 ml	Red top	2 hours
TESTOSTERONE	S257	ECLIA	Serum	2 ml	Red top	2 hours
THYROSCREEN(FT4,T SH)	P133	ECLIA	Serum	2 ml	Red top	2 hours
THYROID PANEL	S255	ECLIA	Serum	2ml	Red top	2 hours
TSH	S259	ECLIA	Serum	2ml	Red top	2 hours
T3	S254	ECLIA	Serum	2ml	Red top	2 Hours
T4	S253	ECLIA	Serum	2ml	Red top	2 Hours
CORTISOL	S149	ECLIA	Serum	2ml	Red top	2 hours
PCT	S220	ELFA	SERUM	2ML	Red top	2hrs
FERRITIN	S166	ELFA	SERUM	2ml	Red Top	2 Hours
VITAMIN D25 HYDROXY	S296	ECLIA	Serum	2ml	Red top	2 Hours

## MICROBIOLOGY

<u>TEST NAME</u>	<u>TEST CODE</u>	<u>METHOD</u>	<u>SPECIMEN</u>	<u>QTY</u>	<u>SAMPLE COLLECTION</u>	<u>REPORT</u>
FUNGAL CULTURE	S175	MANUAL	FLUID	2-5ML	Sterile container	48HRS
FUNGAL CUTURE BLOOD	S176	MANUAL	FLUID	2-5ML	Blood culture bottle(available at lab)	48HRS
FUNGAL SMEAR(KOH)	S177	MICROSCOPY	FLUID,URINE,SPUTUM	5ML	Sterile container	2HRS
CERVICAL SWAB CULTURE	S144	MANUAL	CERVICAL SWAB	-	Sterile container	48HRS
CATHETER TIP CULTURE	S140	MANUAL	CATHETER TIP	-	Sterile container	48HRS
BLOOD CULTURE ( PYOGENICS)	S128	MANUAL	BLOOD CULTURE BOTTLE	2-5ML	Blood culture bottle(available at lab)	7D
BLOOD CULTURE ( ENTERIC )	S126	MANUAL	BLOOD CULTURE BOTTLE	2-5ML	Blood culture bottle(available at lab)	7D
BLOOD CULTURE (PYOGENIC AND ENTERIC )	S127	MANUAL	BLOOD CULTURE BOTTLE	2-5ML	Sterile container	7D
BIOLOGICAL INDICATOR CULTURE	S125	MANUAL	BAL FLUID	-	Sterile container	48HRS
BODY FLUID FOR CULTURE	S115	MANUAL	BODY FLUID	2-5ML	Sterile container	48HRS
BAL FOR CULTURE	S118	MANUAL	BAL FLUID	2-5ML	Sterile container	48HRS
CSF/BODY FLUID , BaL CULTURE (BACT ALERT )	S114	BACT ALERT	BODY FLUID	5-10ML	CULTURE BOTTLES	48HRS
AIR CULTURE	S106	MANUAL	CULTURE PLATE	-	CULTURE PLATE	24HRS
Acid fast bacilli smear(ZNstain)	S302	ZN	Sputum/pus/fluidETC	2-5 ml	Sterile container	Same day
Bactec Blood culture	S129	bactec	Blood inoculated in Bactec bottles	10 ml	Bactec bottle (Pead/ Adult)	48 hrs & 5th day

blood culture & sensitivity Manual	S129	manual	Blood	7-10 ml	Blood culture bottle	48 hrs & 7 th day
CSF C/S BACT ALERT	S155	BACTEC	CSF	5 ML	STERILE CONTAINER	2 DAYS
SPUTUM CULTURE	S242	Manual	SPUTUM	5 ML	STERILE CONTAINER	2 DAYS
STOOL C/S	S243	Manual	STOOL	5 g	STERILE CONTAINER	2 DAYS
URINE C/S	S277	Manual	URINE	5 ML	STERILE CONTAINER	2 DAYS
GRAM STAIN	S182	MANUAL	SPUTUM/PUS /FLUID	5 ML	STERILE CONTAINER	SAME DAY
HIC (S/C) SURVEILLANCE	S158/H100	MICROSCOPY	10 SWABS	SWABS		2 DAYS
THROAT SWAB C/S	S258	Manual	THROAT SWAB	THROAT SWAB	STERILE CONTAINER	2 DAYS
Drinking water culture	S298	manual	Drinking water	150 ml	Sterile container	3 days
India ink	S194	Microscopy	CSF	2 ML	STERILE CONTAINER	SAME DAY
OT SURVEILLANCE (10 SWAB +1 AC )	S193	manual	Swabs , air sample	-	swabs	48 hrs
VAGINAL SWAB FOR C/S	S293	MANUAL	VAGINAL SWAB	-	Sterile container	48HRS
Umbilical line tip c/s	S269	MANUAL	Umbilical line tip	-	Sterile container	48HRS
Tracheal aspirate for c/s	S263	MANUAL	Tracheal aspirate	-	Sterile container	48HRS
Swab c/s	S252	MANUAL	swab	-	Sterile container	48HRS
Suction tip c/s	S251	MANUAL	Suction tip	-	Sterile container	48HRS
Semen c/s	S234	MANUAL	semen	-	Sterile container	48HRS
Nasal secretion for c/s	S211	MANUAL	Nasal secretion	-	Sterile container	48HRS
KOH PREPARATION	S197	MICROSCOPY	FLUID, URINE ,SPUTUM ETC	2ML	Sterile container	2HRS
GASTRIC ASPIRATE C/S	S180	MANUAL	GASTRIC ASPIRATE	-	Sterile container	48HRS

## SEROLOGY

<u>TEST NAME</u>	<u>TEST CODE</u>	<u>METHOD</u>	<u>SPECIM EN</u>	<u>QTY</u>	<u>SAMPLE COLLECTION</u>	<u>REPORT</u>
Antistreptolysin (ASO)	S112	latex agglutination	SERUM	2 ML	RED TOP	2 HOURS
C Reactive Protein (quantitative)	S154	Turbimetry	SERUM	2 ML	RED TOP	2 HOURS
Dengue NS1 antigen+Antibody	S159	Rapid	SERUM	2 ML	RED TOP	2 HOURS
HBsAg	S182	Rapid	SERUM	2 ML	RED TOP	2 HOURS
HCV(RAPID)	S307	Rapid	SERUM	2 ML	RED TOP	2 HOURS
HIV(rapid)	S192	Rapid	SERUM	2 ML	RED TOP+ Consent form	2 HOURS
HAV IgM	S184	Rapid	SERUM	2 ML	RED TOP	2 HOURS
HEV IgM	S189	Rapid	SERUM	2 ML	RED TOP	2 HOURS
Malarial Antigen	S206	Rapid	EDTA Whole Blood	2 ML	EDTA	2 HOURS
URINE PREGNANCY TEST(RAPID)	S290	Rapid	URINE	2 ML	STERILE CONTAINER	1HOURS
Typhi IgM & IgG	S266	Rapid	SERUM	2ML	RED TOP	2 HOURS
Rheumatoid Factor (quantitative)	S227	Turbidimetry	SERUM	2 ML	RED TOP	2 HOURS
Rheumatoid Factor (qualitative)	S226	LATEX SLIDE AGGLUTINATION	SERUM	2 ML	RED TOP	2 HOURS
Skin test- Mantoux	S208	Manual	N/A	N/A	N/A	3 <sup>rd</sup> day
VDRL	S294	RAPID	SERUM	2 ML	RED TOP	2 Hours
VDRL Titre	S295	Flocculation reaction	Serum	2 ml	Red top	2 Hours
Widal Test	S299	Agglutination test	serum	2 ml	Red Top	2 hours